If you detect any defect in the product please contact your dealer or the manufacturer. Should there are particular problems which are not covered sufficiently by these instructions, the matter should be referred to the manufacturer.

Use only milling systems that are recommended by their manufacturers for the processing of zirconium oxide. Milling systems need to be well calibrated for best results. All systems are not alike and can produce adverse results if the minimum thickness is not followed. Because the material shrinks during sintering, it is vital to take into account the appropriate shrinkage factor during milling to ensure the precise fit of the restoration. Each blank is labeled with the specific shrinkage factor that has to be applied.

WARNING: Dust from the milling or grinding process or from pre-sintering manual adjustments may cause irritation of eyes, mucous membranes, skin or damage to lungs. Processing may therefore be performed only with a properly functioning aspirator, safety glasses and an approved dust mask.

The manufacturer accepts no liability for defective results since it has no influence on the milling process.
NOT FOR USE IN THE UNITED STATES ONLY

SINTERING
The sintering furnaces made of priti[®]multidisc ZrO₂ Opaque, priti[®]multidisc ZrO₂ Translucent and priti[®]multidisc ZrO₂ Extra Translucent must go through a final sintering cycle prior to finishing. The sinter firing must only be conducted in high-temperature ovens authorized for this purpose. You may use the firing process specified by pridentia[®]. Please also note the oven manufacturer's specifications. All standard programmable dental laboratory ovens may be used. Sintering ovens vary in their performance, so it is highly recommended the users calibrate periodically to ensure the recommended cycle is carried out correctly. Follow the manufacturer's recommended calibration instructions.

WARNING: Sintering furnaces must be located in a fireproof well-ventilated area. Do not open the oven or remove the sintered restoration before the furnace has cooled down sufficiently for safe handling to avoid the risk of burns.

Recommended firing program:
From room temperature to 1.450 °C at 10 °C/min, hold time of 2 hours at firing temperature, allow to cool to room temperature at 10 °C/min. Natural cooling begins from approx. 600 °C depending on the oven used.
Rapid sintering (optional):
Rapid sintering may be conducted subject to the following conditions. Crowns and bridges with up to 3 units: From room temperature to 1.500 °C at 10 °C/min, hold time of 30 minutes, allow to cool to room temperature at 40 °C/min.

WARNING: Always comply with the above sinter cycles, otherwise the material may become weakened and may break in the patient's oral cavity.
You are strongly advised not to use coloring liquids as this can have a negative effect on translucency or shade.

If restorations have been fabricated in the sintering oven using liquid stains, a cleaning firing with zirconium oxide powder is recommended before priti[®]multidisc ZrO₂ Opaque, priti[®]multidisc ZrO₂ Translucent and priti[®]multidisc ZrO₂ Extra Translucent is sintered in the same oven.

CAUTION: To avoid undesirable color change when sintering pre-shaded zirconium oxide, it is strongly advised to insert a spacer (made of zirconia) of at least 1 mm height between the lid and the tray to ensure that the air can circulate.

Check the sintered product for discoloration, fissures and cracks after processing.
WARNING: Processing a cracked, fractured, broken or discolored product may result in a defective restoration with a risk of injury for patients.

If you detect any defect in the product please contact your dealer or the manufacturer. Should there are particular problems which are not covered sufficiently by these instructions, these matters should be referred to the manufacturer.
NOT FOR USE IN THE UNITED STATES ONLY

INTRODUCTION
Please read these instructions completely and thoroughly before using the product and observe them at all times. Improper use of the product and failure to observe the information provided may affect quality and reduce the lifetime of the dental prostheses. The product must be used in accordance with these instructions and only for the specified purpose. Where it is designed, pridentia[®] GmbH shall not be liable for any consequential damages or injury arising from any use or misuse of this product. In using pridentia[®] GmbH's product you assume the responsibilities of ownership and operation, and you agree to hold pridentia[®] GmbH harmless for any incident or action involving the use of pridentia[®] GmbH's product.

Observe your safe location for further reference over the entire lifetime of the product. Pass on the information provided to any future owner, processor or user of this product or further processed products made from it.

WARNING: Always comply with the above sinter cycles, otherwise the material may become weakened and may break in the patient's oral cavity.
You are strongly advised not to use coloring liquids as this can have a negative effect on translucency or shade.

If restorations have been fabricated in the sintering oven using liquid stains, a cleaning firing with zirconium oxide powder is recommended before priti[®]multidisc ZrO₂ Opaque, priti[®]multidisc ZrO₂ Translucent and priti[®]multidisc ZrO₂ Extra Translucent is sintered in the same oven.

CAUTION: To avoid undesirable color change when sintering pre-shaded zirconium oxide, it is strongly advised to insert a spacer (made of zirconia) of at least 1 mm height between the lid and the tray to ensure that the air can circulate.

Check the sintered product for discoloration, fissures and cracks after processing.
WARNING: Processing a cracked, fractured, broken or discolored product may result in a defective restoration with a risk of injury for patients.

If you detect any defect in the product please contact your dealer or the manufacturer. Should there are particular problems which are not covered sufficiently by these instructions, these matters should be referred to the manufacturer.
NOT FOR USE IN THE UNITED STATES ONLY

ADJUSTMENTS
If adjustments to the densely-sintered restorations are required, these must only be conducted using water-cooled diamond-tipped grinding tools in order to avoid damage of the material by overheating or excessive force on the surface of the restoration. Never use milling tools as they will damage the surface of the restoration.

Basic rules for handling sintered materials:

- work at low pressure,
- use diamond-tipped grinding tools in good condition only,
- avoid sharp edges at the finished product to prevent patient injuries,
- do not machine interdental connectors,
- prior to cementing the restoration in place, clean it and sandblast the internal surfaces in accordance with the following parameters: jet pressure 1 bar, particle size < 50 µm and distance approx. 10 mm.

CERAMIC VENEERS
All veneer ceramics recommended for zirconium oxide ceramics may be used.

FIXING
priti[®]multidisc ZrO₂ Opaque, priti[®]multidisc ZrO₂ Translucent and priti[®]multidisc ZrO₂ Extra Translucent restorations may be cemented conventionally with a zinc phosphate or glass ionomer cement or with self-adhesive fixing and full-adhesive material. Ensure that surfaces are retentive and that stump height is at least 3 mm.

WARNING: Do not use a final restoration if there are any cracked, fractured or broken areas or discolorations visible on the product. If any defect is identified, the framework must not be used as a patient. Observe your safe location for further reference over the entire lifetime of the product. Pass on the information provided to any future owner, processor or user of this product or further processed products made from it.

CONTRAINDICATIONS
WARNING: Do not use the restoration

- for the fabrication of implants
- for the fabrication of abutments
- on patients with parafunctional habits (e.g. bruxism)
- on patients with a known intolerance to components
- with inadequate preparation
- with insufficient spaces available in the patient's mouth
- on patients with inadequate oral hygiene

Otherwise you will risk severe patients' injuries.

CAUTION: Do not store in a damp area. Dampness may affect the product. Do not store near potential contamination sources as they may affect the product.

MATERIAL PROPERTIES
priti[®]multidisc ZrO₂ Opaque, priti[®]multidisc ZrO₂ Translucent and priti[®]multidisc ZrO₂ Extra Translucent round blanks meet the following specifications at completion of the high-density sintering process.

Material characteristics:
Flexural strength: > 1150 MPa
CTE: 10·10⁻⁶ K⁻¹

priti[®]multidisc ZrO₂ Opaque and priti[®]multidisc ZrO₂ Translucent Chemical composition:
Zirconium dioxide (ZrO₂·H₂O): 94.1 – 94.65%
Yttrium oxide (Y₂O₃): 4.65 – 5.95%
Aluminum oxide (Al₂O₃): < 0.4%
Other oxide (Y₂O₃): < 0.7%

priti[®]multidisc ZrO₂ Extra Translucent Chemical composition:
Zirconium dioxide (ZrO₂·H₂O): 92.1 – 92.65%
Yttrium oxide (Y₂O₃): 6.65 – 7.95%
Aluminum oxide (Al₂O₃): < 0.4%
Other oxide (Y₂O₃): < 0.7%

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The proportions of the single components in the total quantity of the components in a blank may vary within the ranges described above; however, the total quantity of components in each single blank will not exceed 100%.

MACHING
As priti[®]multidisc ZrO₂ Opaque, priti[®]multidisc ZrO₂ Translucent and priti[®]multidisc ZrO₂ Extra Translucent is made of a sensitive high performance material always handle it with care. Avoid handling with wet hands. This product may only be managed by trained technicians. The safety warnings in these instructions for use must be complied with. Users bear sole responsibility for the use of the product. The manufacturer accepts no liability for defective results since it has no influence on the machining process.

CONSTRUCTION
Systems must be constructed of the following specifications:

Minimum framework thickness	anterior	0.4 mm	Maryland-bridges	0.4 mm	3-unit bridges	0.6 mm
	posterior	0.6 mm	0.6 mm	0.6 mm		
Connectors	anterior	-	6 mm ²	6 mm ²		
	posterior	-	9 mm ²	9 mm ²		

Framework design: Anatomically shaped (supporting veneer ceramics); fully anatomical

MILLING
CNC (Computer Numerical Control) milling machines, machine tools and processing parameters must be used for milling the blanks. You are strongly advised not to use coolant during the milling process as this may lead to color changes and/or loss of transparency. Check the milled product for discoloration, fissures and cracks after processing.

WARNING: Processing a cracked, fractured, broken or discolored product may result in a defective restoration with a risk of injury for patients.

If you detect any defect in the product please contact your dealer or the manufacturer. Should there are particular problems which are not covered sufficiently by these instructions, the matter should be referred to the manufacturer.
NOT FOR USE IN THE UNITED STATES ONLY

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prifidenta[®] GmbH